

## EBT — iSober 90

Screening operation, printed records, confirmation and escalation

### BEST PRACTICE • DAY-OF-COLLECTION CHECKLIST

After countersigning the consent form, follow the interactive checklist step by step. It is the best way to ensure the correct process and protocols are completed.

### PRODUCT AT A GLANCE

## iSOBER90 breath alcohol device

Use the device with its attached printer and a new, individually packaged mouthpiece for each test.

Photo: Sentechna Canada



### FIELD USE

This card supports trained collectors. It does not replace the iSober 90 manual, current ScreenClear procedure, or eATF.

### BEFORE POWER-ON

- ✓ Verify photo ID and confirm the donor's mouth is clear.
- ✓ Check alkaline AAA power, printer paper, serial number, and calibration record.
- ✓ Confirm that no CAL warning is displayed.

# Power On & Check

Wait for the readiness display before accepting a sample

## FIELD USE

Do not proceed until the device completes its self-check. A device warning is a service issue—not a donor issue.

1

### Power on and wait

Press and hold power. Wait for **bLo**—the internal self-test/readiness display—before accepting a sample.

## READY DISPLAY

**bLo**

Do not collect a sample before the device is ready.



### REMOVE FROM SERVICE

If CAL is displayed, power is low, the device is damaged or dropped, liquid entered the device, or calibration is overdue: remove it from service and contact ScreenClear.

## Obtain & Print

New mouthpiece • steady blow • retain the ticket

### FIELD USE

Use a new mouthpiece for each attempt. The printed ticket is part of the collection record—even for a negative screening result.

**2**

### Fit a new mouthpiece

Open a new sealed, single-use mouthpiece in the donor's presence and attach it firmly. Never reuse a mouthpiece.

**3**

### Obtain the sample

Instruct: "Blow steadily and continuously until you hear two clicks." A steady blow is about 4 seconds; two clicks mean the sample was captured.

### SAMPLE CUE

Steady, continuous blow.

**2 clicks****4**

### Print, photograph and retain

Press print. Take a clear photo of every printed result. Each image uploads automatically and is sent to the required destination. Place each original printed ticket in the provided ScreenClear document folder.

## After the Screening

The screening value decides the next step—not the final outcome

### FIELD USE

Keep the screening ticket with the record. A value at or above 0.020 requires an observed 15-minute wait and confirmation test.

### SCREEN BELOW 0.020

**When:** The printed screening result is below 0.020.

**Action:** Document in the ScreenClear portal. The photo-upload and document-folder requirements apply to every printed result.

### SCREEN 0.020 OR MORE

**When:** The printed screening result is at least 0.020.

**Action:** Record the screening result and start time on the eATF. Begin a **15-minute continuously observed wait**.



### ERR IS NOT A REFUSAL

Repeat the blow. After a second ERR, replace the mouthpiece and repeat. If ERR persists after 3 attempts, stop and contact ScreenClear before making any determination.

# The 15-Minute Wait

A screening result of 0.020 or greater triggers confirmation

## FIELD USE

The confirmation result—not the screening result—controls. Keep the donor continuously observed during the 15-minute wait.

**1**

### Start the observed wait

Record the screening result and start time on the eATF. Start a **15-minute continuously observed wait**.

**2**

### Control the wait

No food, drink, gum, tobacco or vaping, mouthwash, or oral sprays. Restart the timer after belching or regurgitation.

**3**

### Confirm

After 15 minutes, use a new mouthpiece. Perform the confirmation test and print the confirmation ticket.

## Interrupted Wait?

Protect the observation record and obtain direction before continuing

### FIELD USE

Remain factual. Do not attempt to compel the donor or independently determine a refusal.



### LEAVING DURING THE WAIT

Document the facts and contact ScreenClear immediately. Do not attempt to compel the donor or independently make a refusal determination.



Record the time, what occurred, and the point reached in the wait.



Keep completed tickets and the eATF record secure.



Contact ScreenClear before restarting, closing, or changing the process.



### STOP • SECURE • DOCUMENT • CONTACT SCREENCLEAR

1-844-925-3274 | [info@screenclear.ca](mailto:info@screenclear.ca)

Never contact the employer or DER directly.

## Confirmation Result

Document and contact ScreenClear; never contact the employer or DER

### FIELD USE

Attach the printed confirmation ticket to the eATF where required. The collector does not independently direct return to duty.

### < 0.020 — NEGATIVE

**When:** Confirmation below 0.020.

**Action:** Document in the portal. No eATF is required.

### 0.020–0.039 — CAUTION

**When:** Confirmation is at least 0.020 and below 0.040.

**Action:** Complete and submit the eATF. Follow the printed-result photo-upload and document-folder requirements, then contact ScreenClear immediately.

### ≥ 0.040 — ESCALATE

**When:** Confirmation is 0.040 or greater.

**Action:** Complete and submit the eATF. Follow the printed-result photo-upload and document-folder requirements, then contact ScreenClear immediately.

## Before You Close the Record

Keep every ticket, maintain the role boundary, and secure the device

### FIELD USE

Document the collection, retain tickets, secure the device, and contact ScreenClear when direction is needed.



### ROLE BOUNDARY

Do not communicate results to the employer or DER and do not independently direct return to duty. ScreenClear manages employer communication.



Confirm every printed result was photographed, uploaded automatically, and its original ticket is in the provided ScreenClear document folder.



Store the device in its protective case at 5–40°C.



Contact ScreenClear before closing the record if a device or collection concern remains unresolved.